



Ever-Trac Traction Unit

The Ever-Trac Traction Unit is extremely easy to operate and offers all the features expected for lumbar and cervical traction therapy. Supplied with a mains lead, spreader, patient safety release switch, instruction manual and 4 x attachment key wheels.

Features:

- Static, Intermittent, Progressive/Regressive and Cyclic Traction
- Variable speed control
- Minimum and maximum force settings
2.5kgs/5.5lb - 91kgs/200lbs
- Patient safety switch with audible alarm

Specifications:

- Power Supply (Voltage): 100/120/220/240 volts 50/60Hz
Traction Force: 91kg / 200lb maximum
2.5kg/5.5lb Traction force
Compensation: Computerised automatic compensation
Force Display: Actual force/pre-set force - automatic switching display
Traction Speed: Stepless, continuous & adjustable
Therapy Mode: Continuous, Harmonic, Intermittent, Intermittent progressive
Treatment Time: 1 to 99 minutes +- 2%
Current: 1.2A (110V, 120V)
0.6A (220V, 240V)
Power Consumption: 150VA maximum



Code:	Description:
CLA300	Ever-Trac Traction Unit

WHITELEY
HEALTHCARE

